

**WHO Prequalification of In Vitro Diagnostics Programme
PUBLIC REPORT**

**Product: MP Diagnostics HIV Blot 2.2
Number: PQDx 0198-071-00**

Abstract

MP Diagnostics HIV Blot 2.2 with product codes **11030-018** and **11030-036**, manufactured by MP Biomedicals Asia Pacific Pte. Ltd, **CE marked regulatory version**, was accepted for the WHO list of prequalified in vitro diagnostics and was listed on 4 April 2016.

MP Diagnostics HIV Blot 2.2 kit is a qualitative enzyme immunoassay for the in vitro detection of antibodies to HIV-1 and HIV-2 in human serum or plasma. It is intended for use as a supplemental assay on human serum or plasma specimens found repeatedly reactive using EIA. The separated specific HIV-1 viral antigens incorporate onto the strips via electrophoretic and electrotransblot procedures combined with a specific HIV-2 synthetic peptide on the same strip allow for further delineation of the antibody responses to specific viral proteins. Each strip also includes an internal specimen addition control to minimize the risk of false negatives due to operations errors and to ensure the addition of specimens.

The nitrocellulose strips are incorporated with separated bound antigenic proteins from partially unpurified inactivated HIV-1 using electrophoretic blotting, plus a specific HIV-2 synthetic peptide on the same strips. Individual nitrocellulose strips are incubated with diluted specimens and controls. Specific antibodies to HIV-1 and HIV-2 if present in the specimens will bind to the HIV-1 proteins and HIV-2 peptide on the strips. The strips are washed to remove unbound materials. Antibodies that bind specifically to HIV proteins can be visualized using a series of reactions with goat anti-human IgG conjugated with alkaline phosphatase and the substrate BCIP/NBT. This method has the sensitivity to detect marginal amounts of HIV specific antibodies in serum or plasma.

The test kit contains:

Component	18 tests (product code 11030-018)	36 tests (product code 11030-036)
Nitrocellulose Strips: Incorporated with HIV-1 viral lysate, a specific HIV-2 envelope peptide and specimen addition control band.	18 strips	36 strips
Non-reactive control:	1 vial (80µl)	1 vial (80µl)

Inactivated normal human serum non-reactive for hepatitis B surface antigen (HBsAg), antibodies to HIV-1/2, and anti-HCV. Contains sodium azide and thimerosal as preservatives.		
Strong reactive control: Inactivated normal human serum with high titered antibodies to HIV-1 and HIV-2 and non-reactive for HBsAg, and anti-HCV. Contains sodium azide and thimerosal as preservatives.	1 vial (80µl)	1 vial (80µl)
Weak reactive control: Inactivated normal human serum with high titered antibodies to HIV-1 ONLY and non-reactive for HBsAg, anti-HIV-2 and anti-HCV. Contains sodium azide and thimerosal as preservatives.	1 vial (80µl)	1 vial (80µl)
Stock Buffer concentrate (10x): Tris buffer with heat inactivated normal goat serum. Contains thimerosal as presentative.	1 bottle (20ml)	1 bottle (20ml)
Wash buffer concentrate (20x): Tris buffer with Tween-20. Contains thimerosal as preservative.	1 bottle (70ml)	1 bottle (70ml)
Conjugate: Goat anti-human IgG conjugated with alkaline phosphate. Contains thimerosal as preservative.	1 vial (160µl)	1 vial (160µl)
Substrate: Solution of 5-bromo-4-chloro-3-indoyl-phosphate (BCIP) and nitroblue tetrazolium (NBT).	1 bottle (100ml)	1 bottle (100ml)
Blotting powder: Non-fat dry milk	10 packets (1g each)	10 packets (1g each)
Incubation tray: 9 wells each	2 trays	4 trays
Instruction for use:	1 copy	1 copy
Forceps:	1 pair	1 pair

Storage:

The test kit should be stored at 2-8 °C.

Shelf-life:

24 months.

WHO special warning:

WHO reviewed the instructions for use that were current at the time of WHO prequalification, and a number of changes were suggested. Most but not all changes were made by the manufacturer (outstanding comments relate to general layout and nomenclature).

Summary of prequalification status for MP Diagnostics HIV Blot 2.2

	Initial acceptance	
	Date	Outcome
Status on PQ list	4 April 2016	listed
Dossier assessment	10 October 2014	MR
Inspection status	05 June 2015	MR
Laboratory evaluation	N/A	MR

MR: Meets Requirements

NA: Not Applicable

MP Diagnostics HIV Blot 2.2 was accepted for the WHO list of in vitro prequalified diagnostics on the basis of data submitted and publicly available information.

Background information

MP Biomedicals Asia Pacific Pte. Ltd, submitted an application for prequalification of MP Diagnostics HIV Blot 2.2. Based on the established prioritization criteria, MP Diagnostics HIV Blot 2.2 was given priority for prequalification.

Product dossier assessment

MP Biomedicals Asia Pacific Pte. Ltd, submitted a product dossier for MP Diagnostics HIV Blot 2.2 as per the Instructions for compilation of a product dossier (PQDx_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO in accordance with the internal report on the screening and assessment of a product dossier (PQDx_009 v2). Based on the product dossier screening and assessment findings, a recommendation was made to accept the product dossier for MP Diagnostics HIV Blot 2.2 for prequalification.

The manufacturer committed to amend and submit additional documentation on the following issues:

1. The Manufacturer will investigate the effect of transport conditions on HIV Blot 2.2 shelf life including a higher temperature range (expected completion date Q4 2017).
2. The manufacturer will commence testing all material used in the production of positive and negative test kit controls using state-of-the-art methods (i.e. nucleic acid detection).

Manufacturing site inspection

A comprehensive inspection was performed at the site of manufacture (2 Pioneer Place, Singapore) of MP Diagnostics HIV Blot 2.2 in November 2014 as per the Information for manufacturers on prequalification inspection procedures for the sites of manufacture of diagnostics (PQDx_014 v1). The inspection found that the manufacturer had an acceptable quality management system and good manufacturing practices in place that ensured the consistent manufacture of a product of good quality. The manufacturer's responses to the nonconformities found at the time of the inspection were accepted 14 July 2015.

The manufacturer committed conduct the following studies which will be reviewed at the next re-inspection: drop and shock testing for HIV Blot 2.2.

Laboratory evaluation


The objective of the performance laboratory evaluation is to assess the performance and operational characteristics of commercially available in-vitro diagnostics for the purpose of advising the governments of WHO Member States on these issues. In particular, suitability for use in resource-limited settings will be assessed.

Based on the risk level associated with the use of the supplemental assays, the known general performance of supplemental assays and role of the supplemental assays in patient care in resource-limited settings, it was decided that WHO will not conduct performance evaluations of these assays as part of the prequalification assessment process.


Consequently, laboratory evaluation of MP Diagnostics HIV Blot 2.2 was not conducted.

Labelling


- 1. Labels**
- 2. Instructions for use**




Size: 30(W)x27(H)mm
2015-04-20




Size: 30(W)x27(H)mm
2015-04-20




Size: 30(W)x27(H)mm
2015-04-20




Size: 30(W)x27(H)mm
2015-04-20



Size: 80(W)x30(H)mm
2015-08-13




Size: 100(W)x42(H)mm
2015-08-13



Size: 56(W)x27(H)mm
07.10.09



Size: 100(W)x42(H)mm
2015-04-20



Size: 65(W)x45(H)mm
2015-08-05

HIV Blot 2.2

MP Diagnostics HIV BLOT 2.2

WESTERN BLOT ASSAY

- EN** For detection and identification of IgG antibodies to Human Immunodeficiency Virus Type 1(HIV-1) and Type 2 (HIV-2)
- DE** Für den Nachweis und die Identifizierung von IgG-Antikörpern gegen das Humane Immundefizienz-Virus Typ 1 (HIV-1) und Typ 2 (HIV-2)
- FR** Pour la détection et l'identification des anticorps IgG du virus d'immunodéficience humaine de type 1 (VIH-1) et de type 2 (VIH-2)
- IT** Per la determinazione e l'identificazione degli anticorpi IgG del virus dell'immunodeficienza umana di tipo 1 (HIV-1) e di tipo 2 (HIV-2)
- ES** Para la detección e identificación de anticuerpos IgG del virus de inmunodeficiencia humana tipo 1 (VIH-1) y tipo 2 (VIH-2)
- PT** Para a detecção e identificação de anticorpos IgG contra o vírus da imunodeficiência humana tipo 1 (HIV-1) e tipo 2 (HIV-2)

ANTIGEN STRIPS	
CONTROL -	
CONTROL +	
CONTROL WEAK	
BUF STOCK 10x	
BUF WASH 20x	
CONJUGATE	
SUBS BCIP / NBT	
POWDER	BLOTTING


Σ_{18}	Σ_{36}
18 strips	36 strips
80µl x 1	80µl x 1
80µl x 1	80µl x 1
80µl x 1	80µl x 1
20ml x 1	20ml x 1
70ml x 1	70ml x 1
160µl x 1	160µl x 1
100ml x 1	100ml x 1
1g x 10	1g x 10

MP Diagnostics HIV BLOT 2.2

WESTERN BLOT ASSAY

CE
0123

IVD



EC REP

Medical Technology Promedt Consulting GmbH
Altenhofstrasse 80, D-66386 St. Ingbert Germany
Tel. No.: +49-68 94-58 1020 Fax. No.: +49-68 94-58 1021 Email: info@mt-procons.com

Colour Code:



Size: 280mm(W)x83mm(H) 2015-07-22

LAE 0011-3

ANTIGEN STRIPS

contains 100%
Nitrocellulose



Danger

H228 P210 EUH210
P280

BUF STOCK 10x

contains 0.1%
Thimerosal

BUF WASH 20x



Warning

H373 P260 EUH210
P501

HAZ LAE0011-0

HIV Blot 2.2



P485



Black

size: 165mm(W)x83mm(H)

2015-08-13

QUALITY CONTROL

NOTE: Developed strips must be completely dry to avoid misinterpretation.

The presence or absence of antibodies to HIV-1 sample is determined by comparing each nitrocellulose strip to the assay control strips tested with the NON-REACTIVE, STRONG REACTIVE and WEAK REACTIVE controls.

Figure 1a is suggested as an aid to identify the various bands developed on the STRONG REACTIVE Control strip. The Strong Reactive Control as provided in the kit may contain relatively low titer of anti-p55 and anti-p39; as a result, p55 and p39 band for the Strong Reactive Control may appear faintly on the assayed strips. This has no impact on the performance of HIV Blot 2.2 strips in detecting anti-p55 and anti-p39 present in the specimens, as each lot of strip contains sufficient amount of p55 and p39 antigens.

PLEASE NOTE: The numbered end of the strips should be placed at the bottom as shown in the Figure, i.e. the gp120/gp160 bands are the furthest away from the numbered end.

MOLECULAR WEIGHT	GENE	ANTIGEN	DESCRIPTION
gp 160	ENV	Polymeric form of gp41	Broad diffuse glycoprotein
gp 120	ENV	Outer membrane	Diffuse glycoprotein
p66	POL	Reverse Transcriptase	Discreet band
p55	GAG	Precursor protein	Discreet band
p51	POL	Reverse Transcriptase	Discreet band just below p55
p39	GAG	Fragment of p55	Discreet band
gp41	ENV	Transmembrane	Diffuse glycoprotein
p31	POL	Endonuclease	Doublet
p24	GAG	Core protein	Broad band
p17	GAG	Core protein	Broad band

Some of the different antigens mentioned in the Table above are derived from the same precursor protein and may have overlapping epitopes. This should be considered when interpreting the pattern, for example:-

- It is unlikely to detect gp41 in the absence of gp160 because the gp160 is the polymeric form of gp41 and the concentration of gp160 is higher than gp41 on the MP Diagnostics HIV BLOT 2.2. The gp41 appears as a diffuse band. Any sharp and discreet band at the gp41 region should not be interpreted as gp41 band. Many non-HIV infected and normal specimens are found to be reactive to this non-HIV antigen which is likely to originate from the human cell line used to grow the HIV virus.
- p55 is the precursor for p24 and p17. The p55 band is generally detected when there is strong reactivity to p24 and/or p17, it normally appears as a thin band just above p51 band, sometimes these two bands are indistinguishable and may appear as a single band. The bands seen as p42 and p39 are both GAG fragments and should not be interpreted as gp41 (ENV).

- p24 protein is abundant in HIV Blot 2.2 strip. For seroconverting specimens, it is well established that anti-p24 is the first to appear on Western Blot assays. Appearance of p24 band in HIV infected patients would fulfil the positive interpretation criteria for gag protein by WHO, CDC and other international criteria.
- The POL bands p66, p51 and p31 are generally detected simultaneously. However the sensitivity of p66 and p31 are greater than that of p51.
- HIV-2 cross reactivity is variable but typically shows reactivity with GAG and/or POL antigens. However, there can be cross reactivity with the gp160 band in some cases, but rarely with gp41.
- There is also a high molecular weight band around 160KD that is presumed to be a GAG-POL precursor protein. This is seen with some high titered HIV-2 or indeterminate (GAG Reactive Only) sera but the band pattern is a sharp discreet band which is different from the diffuse band of ENV gp160.
- Appearance of single band near p51/p55 is probably an HLA related reactivity (p56), not specific for HIV-1.
- Appearance of p39 and/or p42 without p24 or p17 should not be interpreted.
- Appearance of p66 alone is not HIV-1 specific, but is most likely a reactivity with the host cell proteins (p68).

The interpretation process involves the following:-

- Validate that the serum control band is visible. If the control is negative, the results should be considered invalid as this indicates a technical error such as not adding sample, conjugate or substrate.
- Identify the molecular weight of each band of the test strip using the STRONG and/or WEAK REACTIVE Control strips as a guide.
- Interpretation of the test strip is then based on the detection of specific band patterns as recommended by the appropriate authorities (i.e. Health Ministry, World Health Organization, etc.)

Specific guidelines for interpretation may differ depending on the local policies. MP Diagnostics recommends following the accepted policy to be in accordance with local regulations.

We recommended the following guidelines for the interpretation of the MP Diagnostics HIV Blot 2.2. Results should be recorded for each band detected, result should be interpreted as NEGATIVE, POSITIVE or INDETERMINATE.

PATTERN	INTERPRETATION
No viral specific bands present	NEGATIVE
Detection of p17 antibodies ONLY, no other bands	NEGATIVE
Detection of 2 ENV (gp160/gp41 and gp120) and 1 GAG (p17, p24, p55) or 1 POL (p31, p51, p66)	HIV-1 POSITIVE
Detection of 2 ENV (gp160/gp41 and gp120) and 1 GAG (p17, p24, p55) or 1 POL (p31, p51, p66) and HIV-2 specific band is visible	HIV-1 POSITIVE with HIV-2 INDICATED

Any viral specific bands present but pattern does not meet criteria for POSITIVE	INDETERMINATE ²
Any viral specific bands present but pattern does not meet criteria for POSITIVE but HIV-2 specific band is visible.	INDETERMINATE ² with HIV-2 INDICATED

***INTERPRETATION OF RESULTS FOR INDETERMINATE**

INDETERMINATE results should not be used as the basis for diagnosis of HIV-1 infection. Based on the fact that most persons with an initial INDETERMINATE result who are infected with HIV-1 will develop detectable HIV antibodies within 1 month, US CDC (2001) recommended such persons be re-tested for HIV-1 infection \geq 1 month later. Persons with continued INDETERMINATE results after 1 month are unlikely to be HIV-infected unless recent HIV exposure is suspected.

Based on a recent study of Fiebig *et al* (2003), although the window period for Western Blot in the case of a primary HIV-1 infection could be as long as 22 days, the progression from an INDETERMINATE blot to a full POSITIVE profile took no longer than 8 days. In addition, this laboratory stage of having Western Blot INDETERMINATE was always accompanied with detectable RNA of HIV-1 with cases of true infection. Conversely, no seroconversion was evident in follow-up studies of individuals having screened positive and Western Blot INDETERMINATE results, once confirmed as negative by PCR methods (Sethoe *et al*, 1995). Therefore, it is reasonable to consider persons having Western Blot INDETERMINATE results but additionally tested negative by a RNA test as unlikely to be HIV-infected, especially when the tested individuals are known as not having any risk factor associated with exposure.

In particular, persons having Western Blot INDETERMINATE results derived from a test algorithm using fourth generation ELISAs as the primary screen test should additionally be tested for viral RNA using a molecular-base test such as RT-PCR with primer sets covering HIV-1/2/O. If necessary, a follow-up should be conducted with an additional supplemental assay on a second specimen collected 1 month later. The unique design of fourth generation ELISAs is for a simultaneous detection of both antigen and antibody. Consequently specimens identified as positive by a fourth generation ELISA should contain either antibody or antigen or both. Although more than 95% of those cases of true positive identified by a fourth generation ELISA were anti-HIV related and verifiable (confirmed) by Western Blot (Ly *et al.*, 2000), a supplemental test using RT-PCR for viral RNA detection appeared unavoidable for the small portion of reactivity relating to p24 antigen. Again, persons without any risk of exposure are unlikely HIV-infected, if identified as positive by a fourth generation ELISA accompanied by a Western Blot INDETERMINATE but the findings could not be further supported by a POSITIVE result using a RNA test with primer sets covering HIV-1/2/O.

However, nucleic acid tests (NAT) for HIV DNA or RNA were not approved for diagnostic purpose by the relevant authorities (US CDC, 2001; Constantine & Zink, 2005) until very recently. To date, only one RNA qualitative assay has been approved by the US FDA for diagnosis

of primary and acute infection of HIV-1. Therefore, test algorithms recommended by the US CDC (2001) and WHO (2004) are yet to be updated, and NAT are yet to be included as methods for resolving INDETERMINATE Western Blot results. Nevertheless, US CDC (2001) acknowledged that when in consultation with clinical and infection status among persons with an initial INDETERMINATE Western Blot.

LIMITATION OF THE METHOD

Detection of antibodies to HIV-1 does not constitute a diagnosis of Acquired Immune Deficiency Syndrome (AIDS). A NEGATIVE BLOT is not a guarantee that the causative agent for AIDS is not present. Although a HIV-1 positive test result by Western Blot indicates infection with the virus, a diagnosis of AIDS can only be made clinically if a person meets the case definition of AIDS established by the Center for Disease Control (USA), the World Health Organization or other relevant authorities.

It is known that persons who have recently seroconverted may display incomplete pattern but increase reactivity (both number and intensity of bands) occurs when followed for a period of two to six months. Most blots with POSITIVE results will have other viral specific bands present.

INDETERMINATE results should not be used as the basis for diagnosis of HIV-1 infection. It is recommended that all INDETERMINATE blots be repeated using the original specimen and sequential samples. Blood donors with an INDETERMINATE blot should be re-tested using a fresh specimen after one month (US CDC, 2001). In addition, antibodies to p24 and p31 are known to decrease during the course of AIDS leading to a shift in blot interpretation from POSITIVE to INDETERMINATE. Interpretation of results should then be based on subsequent blot testing and clinical evaluations in such situations.

Due to its highly specific nature, NON-REACTIVITY of samples with HIV-2 specific envelope peptide on an Indeterminate viral blot, does not exclude the possibility of infection with other strains of HIV-2.

Samples that are indicated as HIV-2 infections should be further tested with specific HIV-2 supplemental assays.

SPECIFIC PERFORMANCE CHARACTERISTICS

SENSITIVITY

HIV POSITIVE SAMPLES

The sensitivity of MP Diagnostics HIV Blot 2.2 was evaluated using 209 HIV-1 & 108 HIV-2 positive samples which were well characterized and commercially available. The HIV Blot 2.2 performance was evaluated and compared with established Western blots for HIV-1 and HIV-2.

BIBLIOGRAPHY

- V.C.W.Tsang, K. Hancock, M. Wilson. D.F. Palmer, S. Whaley, J.S. Mc Dougal, and S. Kennedy. March 1985. Developmental Procedure : Enzyme-linked Immuno-electro-transfer Blot technique for HTLV-III/LAV antibodies; CDC, Atlanta.
- H. Towbin, T. Staehlin, and J. Gordon. 1979. Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. Proc. Natl. Acad. Sci.. USA 76: 4350-4354.
- J. Schupbach, M. Popovic, R. V. Gilden. M.A. Gonda, M. G. Sargadharan and R. C. Gallo. 1984. Serological Analysis of subgroup of Human T-Lymphotropic retroviruses (HTLV-III) associated with AIDS. Science 224, 503-505.
- M. G. Sargadharan, M. Popovic, L. Bruch, J. Schupbach and R. C. Gallo. 1984. Antibodies reactive with human T-Lymphotropic retroviruses (HTLV-III) in the serum of patients with AIDS. Science 224, 506-608.
- CDC. 1985. "Provisional public health service inter-agency recommendations for screening donated blood and plasma for antibody to the virus causing Acquired Immune Deficiency Syndrome" - United States Morbidity and Mortality Weekly Report 34 (1) :1-5.
- Proposed World Health Organization 1990 criteria for interpreting results from Western blot assays for HIV-1, HIV-2, and HTLV-I/HTLV-II, Weekly Epidemiological Record 65(37), 281-283.
- F. Clavel, D. Guetard., F. Brun-Vezinet, et al. 1986 Isolation of a new human retrovirus from West African patients with AIDS. Science; 233:343-346.
- F. Clavel., 1987. HIV-2, the West African AIDS virus. AIDS 1:135-140.
- R.S. Tedder, A. Hughes, T. Corrah et al 1988. Envelope cross-reactivity in Western Blot for HIV-1 and HIV-2 may not indicate dual infection. Lancet 11:927-930.
- Bottiger B., A. Karlsson, F. Andreasson et al. 1990. Envelope cross-reactivity between Human Immunodeficiency Virus Type 1 and Type 2 detected by different serological methods: Correlation between cross-neutralization and reactivity against the main neutralizing site. J. Virol. 64(7):3492-3499.
- Centers for Disease Control. 2001. Revised Guidelines for HIV Counseling, Testing, and Referral and Revised Recommendations for HIV Screening of Pregnant Women — United States, Morbid. Mortal. Weekly Rep. 50: RR-19.
- Fiebig, E. W., D. J. Wright, B. D. Rawal, P. E. Garrett, R. T. Schumacher, L. Peddada, C. Heldebrandt, R. Smith, A. Conrad, S. H. Kleinman, and M. P. Busch. 2003. Dynamics of HIV viremia and antibody seroconversion in plasma donors: implications for diagnosis and staging of primary HIV infection. AIDS. 17:1871-1879.
- Ly, T. D., C. Edlinger, and A. Vabret. 2000. Contribution of combined detection assays of p24 antigen and anti-human immunodeficiency virus (HIV) antibodies in diagnosis of primary HIV infection by routine testing. J Clin Microbiol. 38:2459-2461.
- Sethoe, S. Y., A. E. Ling, E. H. Sng, E. H. Monteiro, R. K. Chan. 1995. PCR as a confirmatory test for human immunodeficiency virus type 1 infection in individuals with indeterminate western blot (immunoblot) profiles. J Clin Microbiol. 33:3034-3036.
- Constantine, N. T. and H. Zink. 2005. HIV testing technologies after two decades of evolution. Indian J Med Res. 121:519-538.
- World Health Organization. 2004. Guidelines for HIV Diagnosis and monitoring of antiretroviral therapy. Regional Office for South-East Asia, New Delhi, India.

TECHNICAL PROBLEMS / COMPLAINTS

Should there be a technical problem / complaint, please do the following :
 1. Note the kit lot number, the expiry date and the strip lot number.
 2. Retain the kits and the results that were obtained.
 3. Contact the nearest MP Biomedicals office or your local distributor.

Lyme Disease Mixed Titer Performance Panel PTL201 from BBI (n = 10)	HIV-1 negative = 8 HIV-1 positive = 0 Indeterminate = 2	HIV-1 negative = 1 HIV-1 positive = 0 Indeterminate = 9
Tuberculosis samples from BCP (n = 9)	HIV-1 negative = 3 HIV-1 positive = 0 Indeterminate = 6	HIV-1 negative = 1 HIV-1 positive = 0 Indeterminate = 8
Anti-HCV Mixed Titer Performance Panel (PHV 202) samples from BBI (n = 10)	HIV-1 negative = 4 HIV-1 positive = 0 Indeterminate = 6	HIV-1 negative = 0 HIV-1 positive = 0 Indeterminate = 10
SFTS 94, HTLV I/II positive panel from SFTS, France (n = 10)	HIV-1 negative = 9 HIV-1 positive = 0 Indeterminate = 1	HIV-1 negative = 2 HIV-1 positive = 0 Indeterminate = 8
HEV positive sera from Armed Forces Research Institute of Medical Sciences, Thailand (n = 9)	HIV-1 negative = 0 HIV-1 positive = 0 Indeterminate = 9	HIV-1 negative = 0 HIV-1 positive = 0 Indeterminate = 9
H. pylori positive samples from Dr Roost Laboratory (n = 10)	HIV-1 negative = 4 HIV-1 positive = 0 Indeterminate = 6	HIV-1 negative = 1 HIV-1 positive = 0 Indeterminate = 9
Dengue positive samples from Singapore General Hospital (n = 10)	HIV-1 negative = 10 HIV-1 positive = 0 Indeterminate = 0	HIV-1 negative = 10 HIV-1 positive = 0 Indeterminate = 0
Total no.	167	167
No. of true negatives	117	31
No. of false positives	0	0
No. of indeterminates	49	135
No. of true positives	1	1
Specificity	70.48% (117/166); 95% CI (62.92 % - 77.3 %)	18.67% (31/166); 95% CI (13.06 % - 25.45 %)

Table 6: Potential interfering and pregnant women samples (81 samples)

Panel Name/ Source	MP Diagnostics HIV Blot 2.2	HIV-1 Western blot
Rheumatoid Factor samples from BCP Samples with RF values (0-500), (501-999) & (\geq 1000) are 2, 2 & 14, respectively. (n = 18)	HIV-1 negative = 15 HIV-1 positive = 0 Indeterminate = 3	HIV-1 negative = 0 HIV-1 positive = 0 Indeterminate = 18
Lipemic samples from BiosPacific (n = 10)	HIV-1 negative = 5 HIV-1 positive = 0 Indeterminate = 5	HIV-1 negative = 0 HIV-1 positive = 0 Indeterminate = 10

Table 1: HIV-1 positive samples (209 samples)

Panel Name/ Source	MP Diagnostics HIV Blot 2.2	HIV-1 Western blot
HIV Surveillance Panels from BioClinical Partners, Inc., USA (BCP) 5 panels (10 members each) and 1 panel (9 members) with samples from USA, China, Venezuela, Thailand, Cameroon and India. (n = 59)	HIV-1 positive = 58 Indeterminate = 1	HIV-1 positive = 58 Indeterminate = 1
Panel HIV SFTS 94 from Sanguine Nationale Transfusion Societes (SFTS), France 15 members which were HIV-1 Western blot positive as indicated on the SFTS data sheet were used. (n = 15)	HIV-1 positive = 15 Indeterminate = 0	HIV-1 positive = 9 Indeterminate = 6
HIV-1 positive samples from Boston Biomedical Inc., USA (BBI) and Serologicals Inc., USA (n = 38)	HIV-1 positive = 38 Indeterminate = 0	HIV-1 positive = 38 Indeterminate = 0
HIV-1 positive plasma from BBI (n = 50)	HIV-1 positive = 50 Indeterminate = 0	HIV-1 positive = 50 Indeterminate = 0
HIV-1 positive plasma from LifeBiotech AG (n = 47)	HIV-1 positive = 47 Indeterminate = 0	HIV-1 positive = 47 Indeterminate = 0
Total no.	209	209
No. of true positives	208	202
No. of false negatives	0	0
No. of indeterminates	1	7
Sensitivity	99.52% (208/209); 95% CI (97.36% - 99.99%)	96.65% (202/209); 95% CI (93.22% - 98.64%)

Table 2: HIV-2 positive samples (107 samples)

Panel Name/ Source	MP Diagnostics HIV Blot 2.2	HIV-1 Western blot
HIV Surveillance Panels from BCP 2 panels (10 members each) with samples from Ghana and Nigeria (n = 20)	HIV-2 indicative = 17 Indeterminate = 3	HIV-2 positive = 12 Indeterminate = 8

- Ming Guan, Frequency, causes and new challenges of indeterminate results in Western Blot Confirmatory Testing for Antibodies to Human Immunodeficiency Virus. Clinical and Vaccine Immunology, June 2007, Vol.14, No.6, p649-659.
- Reed DL, et al. HIV-1 Western Blot Standardization. European Ligand Assay Society Sept 1995; 43-47.

MP Biomedicals Asia Pacific Pte Ltd.
 2 Pioneer Place
 Singapore 627885
 Tel. No. : + 65 6775 0008
 Fax. No. : + 65 6774 6146
 Email : enquiry_ap@mpbio.com

EC REP Medical Technology Promedt Consulting GmbH
 Altenhofstrasse 80
 D-66386 St. Ingbert
 Germany
 Tel. No. : + 49 68 94 58 1020
 Fax No. : + 49 68 94 58 1021
 Email : info@mt-procons.com

Regional Office:

MP Biomedicals SAS
 Parc d'Innovation, BP 50067
 67402 Illkirch Cedex
 France
 Tel No. : (33) 388 67 4607
 Fax No. : (33) 388 67 5420
 Email: custserv.eur@mpbio.com

* U.S. Patent 5,721,095


Panel HIV SFTS 94 from SFTS, France. 7 members which were HIV-2 Western blot positive as indicated on the SFTS data sheet were used. One of the 7 members is a 1/10 diluted sample. (n = 7)	HIV-2 indicative = 6 Indeterminate = 1	HIV-2 positive = 7 Indeterminate = 0
HIV-2 positive samples from Dr. Oliveira Varnier, Laboratory of Human Retrovirology, Genova, Italy (n = 45)	HIV-2 indicative = 44 Indeterminate = 1	HIV-2 positive = 43 Indeterminate = 2
HIV-2 positive samples from BBI, BCP and Serologicals (n = 10)	HIV-2 indicative = 10 Indeterminate = 0	HIV-2 positive = 10 Indeterminate = 0
HIV-2 Performance Panel PRF 201 (15 samples) and PRZ 202 (10 samples) from BBI (n = 25)	HIV-2 indicative = 25 Indeterminate = 0	HIV-2 positive = 21 Indeterminate = 4
Total no.	107	107
No. of true positives	102	93
No. of false negatives	0	0
No. of indeterminates	5	14
Sensitivity	95.33% (102/107); 95% CI (89.43% - 98.47%)	86.92% (93/107); 95% CI (79.02% - 92.66%)

SEROCONVERSION

This study was conducted by a third party institution, using a total of 15 commercial seroconversion panels (SeraCare & Zeptomatrix) which were qualified according to common technical specifications for IVD medical devices (2009/886/EC). The seroconversion sensitivity of MPD HIV Blot 2.2 and Chiron RIBA HIV-1/HIV-2 SIA are comparable and both assays reacted similarly in the same panel follow up samples. See Table 3.

Table 3: Performance of Kit based on Positives and/or Indeterminates detected in seroconversion panels

Performance	Panels	No of Panels
MPD and Chiron have equal detection of Positives and Indeterminates	PRB965, PRB966, PRB968, PRB969, PRB970, ZMC6243, ZMC6245, ZMC6246, ZMC9019, PRB9032	10



a. Strong Reactive Control (Reactive for HIV-1 and HIV-2)
 b. Weak Reactive Control (Reactive for HIV-1 only).
 c. Non-Reactive Control.
 d. A typical HIV-2 seropositive serum.

TROUBLE SHOOTING CHART

