

Care, cleaning and disinfection of respiratory equipment

in sterile services department¹

Equipment used for respiratory therapy (e.g. items that come into contact with mucous membranes of the patient) is considered semi critical; such items should be cleaned and then be disinfected before connecting to other patient.² High-level disinfection of respiratory equipment takes place after cleaning and is typically accomplished by physical methods or chemical germicides.³

STEPS



1. Perform hand hygiene.



2. Don personal protective equipment (PPE) to protect against splashing, spraying or aerosols.




3. Wash the equipment with detergent and rinse with clean water.



4. Disinfect

a. **Physical heat** (for heat resistant equipment): e.g. steam, hot water (more than 121°C).

b. **Chemicals** (for plastic and other parts that can be damaged by heat): e.g. hydrogen peroxide > 0.5%, ethanol 70-90% or 0.1% sodium hypochlorite.*



5. Dry

a. **Physical equipment** (e.g. a washer, pasteurizer or autoclave) often has a drying feature within the machine.


b. For **chemical methods**, let equipment parts air dry on a clean towel or cloth.



6. Store equipment dry in closed packages.

(Before storing equipment remove PPE- wash hands Doff, discard PPE and perform hand hygiene.)

CYCLE FINISH



* Disinfectant solutions require preparation and should be used in well-ventilated areas away from patients.

- Prepare a fresh cloth or disposable wipe soaked in a compatible disinfectant and wipe the device from top to bottom, avoiding contact with electrical connectors.
- 0.1% sodium hypochlorite (1000 ppm) should only be used if device is known to withstand use of chlorine-based agents and no ammonia cleaning agents or acidic body fluids (e.g. urine) are present on the device.
- Do not use different disinfectant formulations during the same disinfection step, this may produce toxic fumes.

If chemical disinfection is used, rinse with sterile or clean water (i.e. water boiled for 5 minutes and cooled), or filtered water (i.e. water passed through a 0.2 µ filter), followed by an alcohol rinse and forced-air drying.

1 N.B. It is recommended for all involved in sterile services to review OpenWHO course on Decontamination and sterilization of medical devices: <https://openwho.org/courses/IPC-DECON-EN> As well as refer to WHO's Decontamination and reprocessing of medical devices for health-care facilities: <https://apps.who.int/iris/handle/10665/250232>)

2 Reference: <https://www.who.int/publications/i/item/infection-prevention-and-control-of-epidemic-and-pandemic-prone-acute-respiratory-infections-in-health-care>

3 Please refer to quality assurance and monitoring sections in WHO Guideline for Decontamination and reprocessing of medical devices for health-care facilities: <https://apps.who.int/iris/handle/10665/250232>



World Health
Organization

HEALTH
EMERGENCIES
programme