

 Pharm-Ed	Document type : Check-list	Field /Sub-field : Quality insurance and risk management / hand hygiene	Date : 15.02.2015
---	--------------------------------------	--	-----------------------------

CHECK LIST of the equipment and products contained in a WHO kit for local production of alcohol-based handrub solution - FORMULATION I

Preparation equipment		
<input type="checkbox"/>	Graduated vessel 5000 mL PP graduation (unit)	2+1*
<input type="checkbox"/>	Graduated vessel 2000 mL PP graduation (unit)	1+1*
<input type="checkbox"/>	Beaker 600 mL	2
<input type="checkbox"/>	Measuring cylinder 100 mL, tall, hexagonal base, PP (unit)	1+1*
<input type="checkbox"/>	Measuring cylinder 250 mL, tall, hexagonal base, PP (unit)	1+1*
<input type="checkbox"/>	Measuring cylinder 500 mL, tall, hexagonal base, PP (unit)	2+1*
<input type="checkbox"/>	Measuring cylinder 1000 mL, tall, hexagonal base, PP (unit)	1+1*
<input type="checkbox"/>	Plastic funnel diam. 12 cm (unit)	1+1*
<input type="checkbox"/>	Plastic funnel diam. 20 cm (unit)	1+1*
<input type="checkbox"/>	Industrial bucket of 15 litres PP	1
<input type="checkbox"/>	Rectangular water container PE-LD 20 litres with tap PP (unit)	2+1*
<input type="checkbox"/>	Pair of gloves NITRILE non sterile S/M/L (1 box of each, of 150 units)	3
<input type="checkbox"/>	White bottle 100 mL (unit) without logo Ref. M. R. Senser	3000
<input type="checkbox"/>	Screw tops, leak-proof , blue (unit)	3000
<input type="checkbox"/>	Self-adhesive labels format 5cmx5,6cm: printed by B. Braun	6000
Laboratory equipment		
<input type="checkbox"/>	Pair of plastic security glasses	1+1*
<input type="checkbox"/>	Alcoholmeter Gay LUSSAC (scale v/v)	1+1*
<input type="checkbox"/>	Volumetric flask 10 mL (unit)	4+1*
<input type="checkbox"/>	Graduated Pipette 10 mL PP (unit)	1+1*
<input type="checkbox"/>	Graduated Pipette 1 mL PP (unit)	4+1*
<input type="checkbox"/>	Pipettor (unit)	1+1*
<input type="checkbox"/>	Small funnel	1+1*
<input type="checkbox"/>	Quantofix Peroxyde 1000 strips (1x100)	2
Chemical products		
<input type="checkbox"/>	Glycerol 99,5% Ph. Eur. (Container/5 kG)	2
<input type="checkbox"/>	Hydrogen peroxide 3% (Container/5 litres)	6

*= spare equipment in case of breakage